



RamSoft

# PowerServer

**IHE Integration Statement**

# RamSoft

<b>IHE Integration Statement</b>			
<b>Vendor</b>	<b>Product Name</b>	<b>Version</b>	<b>Release Date</b>
RamSoft Inc	PowerServer, Gateway	4.3	April 19, 2007
This product implements all of the transactions specified in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			
<b>Integration Profiles Implemented</b>	<b>Actors Implemented</b>	<b>Options Implemented</b>	
Scheduled Workflow	Image Display	None	
Consistent Presentation of Images	Image Display, Print Composer	None	
Access to Radiology Information	Image Display, Report Reader, Report Repository, External Report Repository Access	None	
Key Image Note	Image Display	None	
Simple Image and Numeric Report	Report Creator, Report Reader, Report Repository, External Report Repository Access, Enterprise Report Repository	None	
Post-Processing Workflow	Image Display	None	
Reporting Workflow	Image Display	None	
Portable Data for Imaging	Portable Media Creator, Image Display, Report Reader, Print Composer, Portable Media Importer	None	
NM Image	Image Display	None	
<b><u>Internet address for vendor's IHE information:</u></b>			
<a href="http://www.ramsoft.com">http://www.ramsoft.com</a>			
<b>Links to Standards Conformance Statements for the Implementation</b>			

# RamSoft

<b>HL7</b>	<a href="http://www.ramsoft.com">http://www.ramsoft.com</a>	
<b>DICOM</b>	<a href="http://www.ramsoft.com">http://www.ramsoft.com</a>	
<b>Links to general information on IHE</b>		
North America: <a href="http://www.ihe.net">http://www.ihe.net</a>	In Europe: <a href="http://www.ihe-europe.net">http://www.ihe-europe.net</a>	In Japan: <a href="http://www.ihe-j.org">http://www.ihe-j.org</a>