

/HE Integration Statement

Vendor	Product Name	Version	Release Date
RamSoft Inc	PowerServer Lite, Gateway Viewer	6.0	May 27, 2015

This product implements all of the transactions specified in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

Integration Profiles Implemented	Actors Implemented	Options Implemented
Access to Radiology Information (ARI)	Image Display	None
	Image Manager/ Image Archive	None
	Report Reader	None
	Report Repository	None
Consistent Presentation of Images (CPI)	Evidence Creator	None
	Image Display	None
	Image Manager/ Image Archive	None
	Print Composer	None
Consistent Time (CT)	Time Client	None
Evidence Documents (ED)	Evidence Creator	None
	Image Display/ Report Creator	None
	Image Manager/ Image Archive	None

Import Reconciliation Workflow (IRWF)	Image Manager/ Image Archive	None
	Importer	Unscheduled Import
Key Image Note (KIN)	Image Display	None
Nuclear Medicine Image (NMI)	Evidence Creator	None
	Image Display	None
	Image Manager/ Image Archive	None
Patient Information Reconciliation (PIR)	Image Manager/ Image Archive	None
	Report Manager	None
Portable Data for Imaging (POI)	Display	None
	Image Display	None
	Portable Media Creator	None
	Portable Media Importer	None
	Print Composer	None
	Report Reader	None
Radiation Exposure Monitoring (REM)	Image Manager/ Image Archive	None
Reporting Workflow (RWF)	Image Manager/ Image Archive	None
	Report Creator	None
	Report Manager	None

	Report Reader	None
--	---------------	------

Scheduled Workflow (SWF)	Evidence Creator	None
	Image Manager/ Image Archive	None
	Image Display	None
Simple Image and Numeric Report (SINR)	Report Creator	None
	Report Manager	None
	Report Reader	None
	Report Repository	None

Internet address for vendor's IHE information: www.ramsoft.com/resources.php

Links to Standards Conformance Statements for the Implementation

DICOM, HL7 | www.ramsoft.com/resources.php

Links to general information on IHE

www.ihe.net